

Fluid Transfer - Extended Double-Y Elements



WMarchitect
single-use solutions

Single-use fluid transfer assemblies

Features and benefits

- Complete product traceability
- Support of internal validation protocols and validation guides available upon request
- Standard configurations available for reduced turn-around times
- Standard assemblies sold non-irradiated, irradiated and sterile assemblies on request
- Use in combination with Watson-Marlow peristaltic pumps for an end-to-end fluid management solution
- Supply security with Watson-Marlow’s global manufacturing and supply network



Technical specifications

	Fluid Transfer - Extended Double-Y Elements
Certification and compliance	ADCF, ISO10993, USP <788>, USP <85>, USP Class VI
Manufacturing	ISO 14644-1 Class 7, ISO 9001
Sterility	Non-irradiated (irradiated or sterile on request)

	Fluid Transfer - Extended Double-Y Elements
Tubing	Silicone, Thermoplastic elastomer
Connector material	Polycarbonate (PC), Polypropylene, PVDF

Order Information



Part Number	Part Name (Description)	Pump	Pump-head
33-0913-111221	Pumpsil Y-Element : 3.2x2.4mm-1m ends w/open-NON	500 series	505L
33-0913-111222	Pumpsil Y-Element : 3.2x2.4mm -1m ends w/TC-NON	500 series	505L
33-0913-111223	Pumpsil Y-Element : 3.2x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0913-111224	Pumpsil Y-Element : 3.2x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0913-111231	Pumpsil Y-Element : 4.8x2.4mm-1m ends w/open-NON	500 series	505L
33-0913-111232	Pumpsil Y-Element : 4.8x2.4mm-1m ends w/TC-NON	500 series	505L
33-0913-111233	Pumpsil Y-Element : 4.8x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0913-111234	Pumpsil Y-Element : 4.8x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0913-111241	Pumpsil Y-Element : 6.4x2.4mm-1m ends w/open-NON	500 series	505L
33-0913-111242	Pumpsil Y-Element : 6.4x2.4mm-1m ends w/TC-NON	500 series	505L
33-0913-111243	Pumpsil Y-Element : 6.4x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0913-111244	Pumpsil Y-Element : 6.4x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0913-111251	Pumpsil Y-Element : 8.0x2.4mm-1m ends w/open-NON	500 series	505L
33-0913-111252	Pumpsil Y-Element : 8.0x2.4mm-1m ends w/TC-NON	500 series	505L
33-0913-111253	Pumpsil Y-Element : 8.0x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0913-111254	Pumpsil Y-Element : 8.0x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0913-111261	Pumpsil Y-Element : 9.6x2.4mm-1m ends w/open-NON	500 series	505L
33-0913-111262	Pumpsil Y-Element : 9.6x2.4mm-1m ends w/TC-NON	500 series	505L
33-0913-111263	Pumpsil Y-Element : 9.6x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0913-111264	Pumpsil Y-Element : 9.6x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0913-111271	Pumpsil Y-Element : 8x4mm-1m ends w/open-NON	600 series	620L
33-0913-111272	Pumpsil Y-Element : 8x4mm-1m ends w/TC-NON	600 series	620L
33-0913-111273	Pumpsil Y-Element : 8x4mm-1m ends w/MPC-NON	600 series	620L
33-0913-111274	Pumpsil Y-Element : 8x4mm-1m ends w/Aseptiquik-NON	600 series	620L
33-0933-111221	Bioprene Y-Element : 3.2x2.4mm-1m ends w/open-NON	500 series	505L
33-0933-111222	Bioprene Y-Element : 3.2x2.4mm-1m ends w/TC-NON	500 series	505L
33-0933-111223	Bioprene Y-Element : 3.2x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0933-111224	Bioprene Y-Element : 3.2x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0933-111231	Bioprene Y-Element : 4.8x2.4mm-1m ends w/open-NON	500 series	505L
33-0933-111232	Bioprene Y-Element : 4.8x2.4mm-1m ends w/TC-NON	500 series	505L
33-0933-111233	Bioprene Y-Element : 4.8x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0933-111234	Bioprene Y-Element : 4.8x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0933-111241	Bioprene Y-Element : 6.4x2.4mm-1m ends w/open-NON	500 series	505L
33-0933-111242	Bioprene Y-Element : 6.4x2.4mm-1m ends w/TC-NON	500 series	505L
33-0933-111243	Bioprene Y-Element : 6.4x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0933-111244	Bioprene Y-Element : 6.4x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0933-111261	Bioprene Y-Element : 9.6x2.4mm-1m ends w/open-NON	500 series	505L
33-0933-111262	Bioprene Y-Element : 9.6x2.4mm-1m ends w/TC-NON	500 series	505L
33-0933-111263	Bioprene Y-Element : 9.6x2.4mm-1m ends w/MPC-NON	500 series	505L
33-0933-111264	Bioprene Y-Element : 9.6x2.4mm-1m ends w/Aseptiquik-NON	500 series	505L
33-0933-111271	Bioprene Y-Element : 8x4mm-1m ends w/open-NON	600 series	620L
33-0933-111272	Bioprene Y-Element : 8x4mm-1m ends w/TC-NON	600 series	620L
33-0933-111273	Bioprene Y-Element : 8x4mm-1m ends w/MPC-NON	600 series	620L
33-0933-111274	Bioprene Y-Element : 8x4mm-1m ends w/Aseptiquik-NON	600 series	620L

Standard configurations have been designed based on the most used components. Assemblies listed are non-irradiated with 1-meter-long extended ends on each side of the corresponding Double-Y element. Triclamp fittings are uncapped and standard size. Detailed drawings with bill of materials available upon request. For any modifications to length, fitting size/type or irradiation requirement please contact your local Watson-Marlow Representative.

Disclaimer: The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. It is the user's responsibility to ensure product suitability for use within their application. Watson-Marlow and WMArchitect are registered trademarks of Watson-Marlow Limited. BioClamp, BioBarb, FlatBioEndCap, BioEndCap, BioValve, BioTube applicator and Q-Clamp are Trademarks of BioPure Technology Limited. Tri-Clamp is a registered trademark of Alfa Laval Corporate AB.



wmfts.com/global
01 August 2025