

Platinum-cured silicone braided hose

Platinum-cured silicone braided hose

Features and benefits

- Manufactured and packed within an ISO 14644 Class 7 cleanroom
- Post sterilised validation to current USP, EU and FDA standards
- Withstands sterilisation and repeat autoclave
- Available in either 7.6 m (25 ft) and 15.2 m (50 ft) coils



Technical specifications

	Platinum-cured silicone braided hose
Bore size	3.2 - 25.4 mm (0.125 - 1 in)
External diameter	8.8 - 34.5 mm (0.345 - 1.36 in)
Length	7.6 - 15.2 m (25 - 50 Ft)
Certification and compliance	ADCF, ISO 10993, ISO Class 7, Ph. Eur. 3.1.9, USP <788>, USP <85>, USP <87>, USP Class VI
Burst pressure	15.1 - 44.3 bar (220 - 650 psi)
Shelf life	5 years
Gamma stability	Max 50 kGy
Autoclave stability	Yes

Materials of construction

	Platinum-cured silicone braided hose
Inner tube	Platinum-cured silicone
Braid layer 1	Polyester
Outer tube	Platinum-cured silicone

Product codes

Part number	I.D. mm (inch)	O.D. mm (inch)	Length m (ft)	Min. bend radius mm (inch)	*Min. burst pressure at 20C/68F bar (psi)
BPSHP0125-C	3.2 (0.125)	8.8 (0.345)	7.6 (25)	25.4 (1.00)	38.0 (551)
BPSHP0188-C	4.8 (0.187)	10.3 (0.407)	7.6 (25)	31.8 (1.25)	44.3 (650)
BPSHP0250-C	6.4 (0.250)	12.4 (0.490)	7.6 (25)	38.0 (1.5)	43.1 (625)
BPSHP0375-C	9.6 (0.375)	16.7 (0.658)	7.6 (25)	44.4 (1.75)	37.9 (550)
BPSHP0500-C	12.7 (0.500)	20.3 (0.800)	7.6 (25)	50.8 (2.00)	34.5 (500)
BPSHP0625-C	15.9 (0.625)	24.5 (0.965)	7.6 (25)	63.5 (2.5)	28.9 (420)
BPSHP0750-C	19.0 (0.750)	27.9 (1.100)	7.6 (25)	76.2 (3.00)	24.1 (350)
BPSHP1000-C	25.4 (1.00)	34.5 (1.360)	7.6 (25)	101.6 (4.00)	15.1 (220)
BPSHP0125-D	3.2 (0.125)	8.8 (0.345)	15.2 (50)	25.4 (1.00)	38.0 (551)
BPSHP0188-D	4.8 (0.187)	10.3 (0.407)	15.2 (50)	31.8 (1.25)	44.3 (650)
BPSHP0250-D	6.4 (0.250)	12.4 (0.490)	15.2 (50)	38.0 (1.5)	43.1 (625)
BPSHP0375-D	9.6 (0.375)	16.7 (0.658)	15.2 (50)	44.4 (1.75)	37.9 (550)
BPSHP0500-D	12.7 (0.500)	20.3 (0.800)	15.2 (50)	50.8 (2.00)	34.5 (500)
BPSHP0625-D	15.9 (0.625)	24.5 (0.965)	15.2 (50)	63.5 (2.5)	28.9 (420)
BPSHP0750-D	19.0 (0.750)	27.9 (1.100)	15.2 (50)	76.2 (3.00)	24.1 (350)
BPSHP1000-D	25.4 (1.00)	34.5 (1.360)	15.2 (50)	101.6 (4.00)	15.1 (220)

*Burst pressure measured in accordance with the newly adopted ISO1402

Disclaimer: The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. It is the user's responsibility to ensure product suitability for use within their application. Watson-Marlow and WMArchitect are registered trademarks of Watson-Marlow Limited. BioClamp, BioBarb, FlatBioEndCap, BioEndCap, BioValve and BioTube applicator are Trademarks of BioPure Technology Limited. Tri-Clamp is a registered trademark of Alfa Laval Corporate AB. Q-Clamp is used under licence.



wmfts.com/global
29 April 2026