

Accusil Elements

Accusil tubing

Features and benefits

- Excellent dispensing accuracy with $\pm 0.5\%$
- 100 % traceability with laser-etched part-number, lot number and use-by date
- Eight different sizes for filling between 0.0034 fl oz (0.1 ml) to more than 8.45 fl oz (250 ml)

Technical specifications

	Accusil Elements
Sterilisation methods	Autoclave, EtO, Gamma, X-ray
Autoclave stability	121 °C 30 minutes
Gamma stability	Max 50kGy
Bore size	0.0625 - 0.315 in (1.6 - 8 mm)
Wall thickness	0.09375 in (2.4 mm)
Colour	Colorless
Transparency	Translucent
Spallation	Low
Operating temperature	-4 °F to 176 °F (-20 °C to 80 °C)
Hardness, shore A (5sec)	58 to 64
Specific gravity	1.15 g/ml
Compression set	$\leq 30\%$
Shelf life	5 years

Materials of construction

	Accusil Elements
Tubing	Platinum-cured silicone
Connector material	PVDF

Product Codes

Bore		Wall		Double segment kynar Y-element
mm	inch	mm	inch	
1.6	1/16	2.4	3/32	84-510-016
3.2	1/8	2.4	3/32	84-510-032
4.8	3/16	2.4	3/32	84-510-048
6.4	1/4	2.4	3/32	84-510-064
8.0	5/16	2.4	3/32	84-510-080
9.6	3/8	2.4	3/32	84-510-096

Disclaimer: The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. It is the users responsibility to ensure product suitability for use within their application. Watson-Marlow and WMArchitect are registered trademarks of Watson-Marlow Limited. Bio Y, BioClamp, BioBarb, FlatBioEndCap, BioEndCap, BioValve and BioTube applicator are Trademarks of BioPure Technology Limited. Tri-Clamp is a registered trademark of Alfa Laval Corporate AB. Q-Clamp is used under license.



wmfts.com/global
28 October 2025