

# Small Silicone RXPX Series Gaskets



*Silicone RXPX series gaskets*

## Features and benefits

- Designed in accordance with ASME-BPE standards
- USP Class VI compliant and animal derived component free (ADCF)
- Compliant with FDA regulations CFR 21 177 2600
- Lot numbering enables product traceability
- ISO 14644-1 Class 7 clean room produced and packed
- Single polyethylene (PE) bag




## Technical specifications

|                              | Small Silicone RXPX Series Gaskets  |
|------------------------------|-------------------------------------|
| Gasket type                  | Unflanged                           |
| Internal diameter ranges     | 10.9 - 16 mm (0.43 - 0.63 in)       |
| External diameter ranges     | 21.6 mm (0.85 in)                   |
| Size                         | 0.5, 0.75 in                        |
| Quantity                     | 25 PCS                              |
| Certification and compliance | ADCF, ISO Class 7, USP Class VI     |
| Operating temperature ranges | -55 °C to 210 °C (-67 °F to 410 °F) |
| Shelf life                   | 5 years                             |
| Gamma stability              | Max 50 kGy                          |

## Materials of construction

|          | Small Silicone RXPX Series Gaskets |
|----------|------------------------------------|
| Material | Platinum-cured silicone            |

Product codes



| Part number   | Size inch | O.D inch (A) | I.D inch (B) | Quantity |
|---------------|-----------|--------------|--------------|----------|
| 42RXPX-050-25 | 1/2       | 0.85 (21.6)  | 0.43 (10.9)  | 25       |
| 42RXPX-075-25 | 3/4       | 0.85 (21.6)  | 0.63 (16.0)  | 25       |

Disclaimer: The information contained in this document is believed to be correct but Watson-Marlow Limited accepts no liability for any errors it contains and reserves the right to alter specifications without notice. It is the user's responsibility to ensure product suitability for use within their application. Watson-Marlow and WMArchitect are registered trademarks of Watson-Marlow Limited. BioClamp, BioBarb, FlatBioEndCap, BioEndCap, BioValve, BioTube applicator and Q-Clamp are Trademarks of BioPure Technology Limited. Tri-Clamp is a registered trademark of Alfa Laval Corporate AB.



wmfts.com/global  
13 October 2025